



891 "Flake"

Independent chambers vacuum mattress



FEATURES

Vacuum mattress with 3 independent chambers for an excellent and fast immobilization. Extremely light and sturdy, 8 handles, restraining belts and a perfect patient encirclement, give to the mattress unique solutions. Unmatchable vacuum durability without any deflection. X-ray transparent mattress.

MAIN FEATURES

| | |
|-----------------------|----------------------|
| Size | Adults |
| Usage | Immobilization |
| Length (cm) | 200 |
| Width (cm) | 90 |
| Weight (kg) | 6,7 |
| Loading capacity (kg) | 150 |
| Trauma | |
| Material | PVC-coated polyester |

SEE ALSO



891/SG

Flake basic independent chambers vacuum mattress



Chesterfield 9114

Vacuum mattress with independent chambers with central welding

ACCESSORIES



888
Aluminium aspiration pump



890-001
Carrying case for vacuum mattress



896/DF
Abs aspiration pump

STANDARD ACCESSORIES



894-001
Aspiration valve for vacuum splints
and mattress
