



891 "Flake"

Independent chambers vacuum mattress



FEATURES

Vacuum mattress with 3 independent chambers for an excellent and fast immobilization. Extremely light and sturdy, 8 handles, restraining belts and a perfect patient encirclement, give to the mattress unique solutions. Unmatchable vacuum durability without any deflection. X-ray transparent mattress.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Length (cm)	200
Width (cm)	90
Weight (kg)	6,7
Loading capacity (kg)	150
Trauma	
Material	PVC-coated polyester

SEE ALSO



891/SG

Flake basic independent chambers vacuum mattress



Chesterfield 9114

Vacuum mattress with independent chambers with central welding

ACCESSORIES



888
Aluminium aspiration pump



890-001
Carrying case for vacuum mattress



896/DF
Abs aspiration pump

STANDARD ACCESSORIES



894-001
Aspiration valve for vacuum splints
and mattress
