



431

Arcos rigid splint set with carrying bag



FEATURES

New rigid splints, made with fire-proof and water proof PVC, allow a correct immobilization of the limbs. They are X-ray translucent. Equipped with 4 Velcro® closures with adjustable length allowing operator to apply them easily and safely. Shaped to strongly and easily fix to ankle and upper part of the leg, they are supplied with foot support integrated.

Available in two sizes: Adult arm/ pediatric leg: lenght cm. 52 - base cm 9

Adult leg: lenght cm 74 - base cm 13 Pediatric limb: lenght 30 cm, base 9 cm

MAIN FEATURES

Trauma	⊌′
Usage	Immobilization

IMAGES



SEE ALSO



425 Arcos long rigid splint for lower limbs



427 Arcos short rigid splint for upper limbs



429 Arcos Splint for pediatric limbs