



427

Arcos short rigid splint for upper limbs



FEATURES

New rigid splints, made with fire-proof and water proof PVC, allow a correct immobilization of the limbs. They are X-ray translucent. Equipped with 4 Velcro® closures with adjustable length allowing operator to apply them easily and safely. Shaped to strongly and easily fix to ankle and upper part of the leg, they are supplied with foot support integrated.

Available in two sizes:

Adult arm/ pediatric leg: length cm. 52 - base cm 9

Adult leg: length cm 74 - base cm 13

MAIN FEATURES

Pediatric	
Trauma	
Length (cm)	52
Width (cm)	24
Weight (kg)	0,7
Usage	Immobilization
Material	PVC

IMAGES



SEE ALSO



430
Arcos rigid splint set with carrying bag



425
Arcos long rigid splint for lower limbs