



## 9124/B


Standard spinal immobilizer, blue colour



### FEATURES

Essential device used to safely immobilize and extricate a patient. Made of vinyl-coated nylon. Complete with a padded pillow, two head straps and carrying bag.

### MAIN FEATURES

Trauma	
Usage	Immobilization
Width (cm)	85
Height (cm)	83
Belts' length (cm)	120
Max circumference	180
Weight (kg)	2,9
Loading capacity (kg)	230
Material	Nylon

---

SEE ALSO

---



**Med Ped 9120**  
Pediatric spinal immobilizer

---