




629 Jelly


Red standard head immobilizer for spinal board



FEATURES

Adjustable to every models of spinal boards. It is made of expanded materials at closed cells covered by a vinyl material.

MAIN FEATURES

Trauma	
Width (cm)	40,5
Depth (cm.)	25,5
Height (cm)	16,5

SEE ALSO



629/SLIM
Red head immobilizer with black
slim base for spinal board

ACCESSORIES



629-007

Slim black plate for art. 629 head immobilizer

ACCESSORY OF




Galaxy 9010

"GALAXY" Yellow spinal board



Galaxy 9012

"GALAXY" Orange spinal board



Galaxy-P 9040

Yellow spinal board with pins



Galaxy-P 9042

Orange spinal board with pins



Nilix 9050

Yellow basic spinal board with pins
