




629 JELLY


Orange standard head immobilizer for spinal board



CARATTERISTICHE

Adjustable to every models of spinal boards it is made of expanded materials at closed cells covered by a vinyl material.

HAUPTEIGENSCHAFTEN

| | |
|-------------|---|
| Trauma |  |
| Breite (cm) | 40,5 |
| Tiefe (cm.) | 25,5 |
| Höhe (cm) | 16,5 |

ZUBEHÖR VON




Galaxy 9010
Spine board, gelb



Galaxy 9012
Spine board, orange



Galaxy 9015
Spine board, blau



Galaxy-P 9040
Tavola spinale gialla con pins



Galaxy-P 9042
Tavola spinale arancio con pins



Nilix 9050
Tavola spinale gialla basic con pins