




## 629 Jelly

---

Orange standard head immobilizer for spinal board




---

### FEATURES

Adjustable to every models of spinal boards. It is made of expanded materials at closed cells covered by a vinyl material.

---

### MAIN FEATURES

Trauma	
Width (cm)	40,5
Depth (cm)	25,5
Height (cm)	16,5

---

## ACCESSORY OF

---




**Galaxy 9010**  
"GALAXY" Yellow spinal board



**Galaxy 9012**  
"GALAXY" Orange spinal board



**Galaxy 9016**  
"GALAXY" Red spinal board



**Galaxy-P 9040**  
Yellow spinal board with pins



**Galaxy-P 9042**  
Orange spinal board with pins



**Nilix 9050**  
Yellow basic spinal board with pins