



# 676 Extra Lite Standard

---

Foldable carrying chair with 4 wheels, footrest and arms



---

## FEATURES

EXTRA LITE Standard is the new Me.Ber. stair chair entirely made of lightweight aluminium alloy frame. The backrest and back carrying handles inclination, along with the unique frontal arms configuration, enable a perfect weights distribution, relieving operators' efforts. While manoeuvring on smooth and plane surfaces, the 200 mm rear and 125 mm front wheels gives exceptional comfort. The structure along with an ergonomic seat and backrest, turnable foot rest with magnetic locks and armrest complete the characteristics of this new rescue device.







Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789.

---

## MAIN FEATURES

---

Type	
Foldable	
Track assisted system	
Height (cm)	101
Width (cm)	55
Length (cm)	78
Max length (cm)	139
Epaisseur pliée (cm)	29
Weight (kg)	13,9
Loading capacity (kg)	230
Folded height (cm)	101
Footboards and arms	
Cardio	
10G Certification	
Material	Aluminium alloy

---

## SEE ALSO

---



**679**  
Narrow foldable sedan chair



**678 Extra Lite Basic**  
Extra Lite Basic, foldable carrying chair

---

## ACCESSORIES

---



**960**

"Extra fixing system" - fixing system for chair Extra



**670-001**

Orange/black head immobilizer for Extra stair chair



**670-002**

Transpire, dust proof black nylon cover for Extra chairs



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**975**

"Ergolift fixing system" - Fixing system for chair Extra Ergolift



**670-004**

Castor with brake diam. 200 mm for Extra chair

---

## STANDARD ACCESSORIES

---



**566**

Black bag with 3 belts in 2 pieces made of black ribbon with orange MeBer logo for Extra chairs

---