



676 Extra Lite Standard

Foldable carrying chair with 4 wheels, footrest and arms







FEATURES

EXTRA LITE Standard is the new Me.Ber. stair chair entirely made of lightweight aluminium alloy frame. The backrest and back carrying handles inclination, along with the unique frontal arms configuration, enable a perfect weights distribution, relieving operators' efforts. While manoeuvring on smooth and plane surfaces, the 200 mm rear and 125 mm front wheels gives exceptional comfort. The structure along with an ergonomic seat and backrest, turnable foot rest with magnetic locks and armrest complete the characteristics of this new rescue device.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789.

MAIN FEATURES

Type	
Foldable	
Track assisted system	
Height (cm)	101
Width (cm)	55
Length (cm)	78
Max length (cm)	139
Epaisseur pliée (cm)	29
Weight (kg)	13,9
Loading capacity (kg)	230
Folded height (cm)	101
Footboards and arms	
Cardio	
10G Certification	
Material	Aluminium alloy

SEE ALSO



679
Narrow foldable sedan chair



678 Extra Lite Basic
Extra Lite Basic, foldable carrying chair

ACCESSORIES



960
"Extra fixing system" - fixing system
for chair Extra



670-001
Orange/black head immobilizer for
Extra stair chair



670-002
Transpire, dust proof black nylon
cover for Extra chairs



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



975
"Ergolift fixing system" - Fixing
system for chair Extra Ergolift



670-004
Castor with brake diam. 200 mm
for Extra chair

STANDARD ACCESSORIES



566
Black bag with 3 belts in 2 pieces
made of black ribbon with orange
MeBer logo for Extra chairs
