



678 Extra Lite Basic

Extra Lite Basic, foldable carrying chair



FEATURES

EXTRA LITE basic is the new Me.Ber. stair chair entirely made of lightweight aluminium alloy frame. The backrest and back carrying handles inclination, along with the unique frontal arms configuration, enable a perfect weights distribution, relieving operators' efforts. While manoeuvring on smooth and plane surfaces, the 200 mm rear and 125 mm front wheels gives exceptional comfort. The structure along with an ergonomic seat and backrest, complete the characteristics of this new rescue device.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

Type	
Foldable	
Track assisted system	
Height (cm)	101
Width (cm)	55
Length (cm)	76
Max length (cm)	139
Epaisseur pliée (cm)	29
Weight (kg)	13,6
Loading capacity (kg)	230
Folded height (cm)	101
Cardio	
10G Certification	
Material	Aluminium alloy

SEE ALSO



679
Narrow foldable sedan chair



Extra Lite Standard 676
Foldable carrying chair with 4 wheels, footrest and arms

ACCESSORIES



670-001
Orange/black head immobilizer for Extra stair chair



960
"Extra fixing system" - fixing system for chair Extra



566
Black bag with 3 belts in 2 pieces made of black ribbon with orange MeBer logo for Extra chairs



975
"Ergolift fixing system" - Fixing system for chair Extra Ergolift



670-004
Castor with brake diam. 200 mm for Extra chair



670-002
Transpire, dust proof black nylon cover for Extra chairs

STANDARD ACCESSORIES



566

Black bag with 3 belts in 2 pieces
made of black ribbon with orange
MeBer logo for Extra chairs
